
1. INTENDED USE

The HemoCue® HbA1c 501 Test Cartridge, together with the HemoCue® HbA1c 501 Analyzer, which are parts of the HemoCue® HbA1c 501 system, provide a convenient method for measuring the percentage of hemoglobin A1c (HbA1c) in blood samples, and anticoagulated venous whole blood samples. The test is for point of care use to monitor glycaemic control in patients with diabetes mellitus. The HemoCue® HbA1c 501 system uses a boronate affinity assay to separate the glycated hemoglobin fraction from the nonglycated fraction. The HemoCue HbA1c 501 System is intended to be used by professionals in laboratories, clinics and hospitals.

2. TEST CARTRIDGE

The HemoCue® HbA1c 501 Test Cartridges are intended to be used in the HemoCue® HbA1c 501 Analyzer only.

Chapter 1 (reagent solution) contains: 20 v/v % boronate affinity bead, surfactant, 0.1 v/v % non-reactive ingredients, buffer.

Chapter 2 (rinsing solution) contains: surfactant, 0.1 v/v % non-reactive ingredients, buffer. Additional items needed for testing: sterile and single-use liquid control, disposable gloves, tissue, cotton balls, gauze.

3. STORAGE AND LIFE USE

Test Cartridge
Store test cartridges in the protective packaging until use. The test cartridge is stable until the expiration date printed on the package when stored at 2 - 32 °C (36° - 90°F). Do not freeze.

Use life
Use the test cartridge directly after opening it (maximum 2 minutes).

Preparation of the Test Cartridge
Remove test cartridge from protective packaging. If refrigerated, allow to reach room temperature for 30 minutes before use. The operating temperature of the test is 17 - 32°C (63 - 90°F).

4. PRECAUTION/WARNING:
• For in vitro diagnostic use only. • If refrigerated, allow sealed pouches reach room temperature 30 minutes before use. • Do not use cartridges after expiry date on pack and box. • Do not use scissors to open the foil pouch. • Handle and dispose of all samples and pipetting following appropriate biohazard procedures. • Wear personal protective equipment when handling all reagents. Samples, quality controls and while operating the HemoCue® HbA1c 501 Analyzer. • Waste material containing patient samples, quality controls or biological products should be considered as hazardous. • Do not touch the optical part of the Cartridge. • Dispose all waste in accordance with applicable national and/or local regulations.

5. PROCEDURE

Please read the relevant Operating Manual for proper use of the system.

Expected values
The American Diabetes Association's (ADA)'s 2014 Clinical Practice Recommendation for diabetes (6) specifies a treatment goal of less than 7 % HbA1c.

Limitation of Procedure:
The HemoCue® HbA1c 501 assay gives accurate and precise results over a total haemoglobin range of 7 to 20 g/dL. Most patients have hemoglobin concentrations within this range. However, patients with severe anemia may have hemoglobin concentrations below 7 g/dL, and patients with polycythemia may have hemoglobin concentrations above 20 g/dL. Patients known to have these conditions should be tested with another method for HbA1c determination.

6. CALIBRATION

HemoCue® HbA1c 501 Analyzer:
The HemoCue® HbA1c 501 Analyzer is factory calibrated. The Analyzer automatically self-adjusts during each warm-up and during each assay. An error message will appear if the analyzer is not able to make appropriate internal adjustments.

HemoCue® HbA1c 501 Test Cartridge:
At Test Cartridges production, each production lot undergoes a thorough analysis and characterization, before being released. The HemoCue® HbA1c 501 system test method is National Glycohemoglobin Standardization Program (NGSP) certified, thus the values of calibration parameters determined to provide an optimal reagent performance, are based on the Diabetes Control and Complications Trial (DCCT) reference method.

Calibration parameters are printed for each of the Test Cartridges on the Test Cartridge label. As soon as the Test Cartridge is inserted into the analyzer, the system automatically recognizes the cartridge code. This accesses the appropriate calibration parameter values (calibration curve) for the particular lot number of reagent packs to use.

7. QUALITY CONTROL
Daily check cartridge and Monthly check cartridge are required to run the HemoCue HbA1c 501 system. The system can be verified according to local guidelines using a commercially available HbA1c control. Follow instructions in HemoCue HbA1c 501 operating manual, see section "System check".

8. IFC STANDARDIZATION

The relationship between HbA1c results from the NGSP network and the IFC network has been evaluated and a master equation has been developed. IFC® results are consistently 1.5– 2 % HbA1c lower throughout the range of values compared to NGSP results. IFC® value = 1.093 x NGSP value + 1.50

As of October 2018, all IFC Network Laboratories are reporting their IFC® HbA1c result in mmol/L or mg/dL. As of October 2018, all HbA1c Network Laboratories are following their IFC® HbA1c result in mmol/L or mg/dL. The equation is expressed as follows: IFC® HbA1c (mmol/mol) = 1.093 x NGSP HbA1c (%) - 23.50

9. PERFORMANCE CHARACTERISTICS
Clinical Correlation
The system accuracy evaluation was performed with 131 fresh venous and capillary blood samples, each sample was

measured by two different analyzers and in duplicate.

The percentage of HbA1c in clinical samples ranging from 4.4 % - 13.9 % HbA1c (venous and capillary) was determined using the HemoCue® HbA1c 501 system and ion exchange high performance liquid chromatography (HPLC) instrument. Results are as follows:

Sample Type	Regression Line	Slope	Intercept	Coefficient of determination (R ²)
Venous	y = 0.9972 x + 0.0495	0.9972	0.0495	0.9850
Capillary	y = 0.9973 x + 0.0297	0.9973	0.0297	0.9890
Total	y = 0.9971 x + 0.0329	0.9971	0.0329	0.9886

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KIRJALLISUUSVIITTEET – The Diabetes Control and Complications Trial Research Group. N. Eng. J. Med. 1993; 329: 977-986. • UK Prospective Diabetes Study Group. Lancet 1998; 352: 837-853. • Baynes, J. W., Burn, H. F., and Goldstein, D. E., et al. National Diabetes Data Group. Report of the Expert Committee on Glycosylated Hemoglobin. Diabetes Care 7 (1984): 602-606. • Koening, R. J., Peterson, C. M., and Kilo, C., et al. Hemoglobin A1c as an indicator of the degree of glucose intolerance in diabetes. Diabetes 25 (1976): 230-232. • Nathan, D. M., Singer, D. E., Hurttal, K., and Goodson, J. D. The clinical information value of the glycosylated hemoglobin assay. N. Eng. J. Med. 310 (1984): 341-346. • American Diabetes Association. Standards of Medical Care in Diabetes - 2012. DIABETES CARE, VOLUME 35, SUPPLEMENT 1, JANUARY 2012. • Larsen, M. L., Horder, M., and Mogensen, E. F. Effect of long-term monitoring of glycosylated hemoglobin levels in insulin-dependent diabetes mellitus. N. Eng. J. Med. 323 (1990): 1021-1025. • Nathan, D. M., Hemoglobin A1c - Infatuation or the real thing? N. Eng. J. Med. 323 (1990): 1062-1063. • HemoCue HbA1c 501 -käyttöohje.

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