Control kit for quality control for the hsCRP parameter on the smart analyser.

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### Kit Size

1 x 1 ml

**Test kit preparation: Allow single test min. 10 minutes to warm to room temperature by placing the test into the test kit rack. Put test kit package back into refrigerator.**

### Summary

hsCRP control kit is intended for use in the clinical laboratory as a control for quantitative determination of hsCRP. It is prepared from a human serum and was found to be non reactive when tested with approved methods for HbsAg, anti HIV1/2 and anti HCV. However this material should be handled as though capable of transmitting infectious disease.

### Storage and Handling

Unopened vials are stable at 2-8°C until the expiration date on the label. Thoroughly mix the content of the vial before each use by gently inverting it. After use please cap vial immediately and put back into refrigerator. Once opened vials are stable until the expiry date when stored tightly capped at 2-8°C.

We recommend not to pipette directly from the original bottle as well as using single-use pipetting tips to avoid contamination. Dispense about 5 drops into a serum cup and pipette out of it. Dispose cup after use.

### Waste Management

Please refer to your local legal requirements.

### Preparation

The controls are liquid and ready to use. For further details on how to use the control kit please refer to the description inside this package insert.

### Control Values

The value is traceable to the IFCC/BCR/CAP reference material for 15 plasma proteins CRM470.

### Warnings and Precautions

- Y Set the sex to "control" and the sample type to "serum" before running the QC serum.
- Y Compare the results with the values indicated in the lot specific value sheet.
- Y If the result is not within the acceptable range repeat measurement. If result is again not within the acceptable range please contact your distributor.

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

## Processing of hsCRP control measurement

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




### ATTENTION!

Do not wipe the capillary or pipetting tip after aspirating the sample! Allow single test min. 10 minutes to warm to room temperature!

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## English

- 1. Preparation of test system**
  - 1.1 Insert RFID card
  - 1.2 Place ERS cuvette in smart test kit rack
  - 1.3 Place ERS cap in smart test kit rack
  - 1.4 Press symbol on touch screen, required information is entered by use of the touch display of the analyser
- 2. Attention: Select serum and control – Measurement on the smart analyser**
  - 2.1 Sex menu: Select control
  - 2.2 Samplytype menu: Select serum
  - 2.3 Instrument ready for start of analysis
- 3. Processing hsCRP control measurement with pipette**
  - 3.1 Dispense some drops control liquid into serum cup. Aspirate 20µl from serum cup.
  - 3.2 Dispense 20µl sample INTO LIQUID in ERS cuvette
  - 3.3 Apply ERS cap firmly onto ERS cuvette
  - 3.4 Place ERS cartridge into smart analyser.
  - 3.5 Start automatic sample processing by pressing the  start button.

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H

M Eurolyser Diagnostica GmbH  
www.eurolyser.com

Technical details subject to change without notice.

Doc-ID: 166c3135-b3d1-4d03-a075-e12244b16461 ; latest revision: DP1.0 ; date of print: 2013-02-25