HemoCue® Glucose 201+

Operating Manual

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
Table of Contents

HemoCue® Glucose 201+ system	5
Components	6
Start up	8
Set-up.....	10
Audio signal, time and date.....	10
Measuring	12
Capillary blood.....	12
Venous or arterial blood	20
Set-up.....	22
QC-test	22
Memory function – Scroll	24
Delete results	26
Printer function	28
Maintenance	30
Daily maintenance	30
Optical parts	32
Troubleshooting Guide.....	34
Specifications.....	46

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GB

Thank you for choosing the HemoCue® Glucose 201+. The system is a fast, simple and reliable method for the determination of glucose in whole blood. The results obtained are comparable to those from standard laboratory methods. Capillary, venous or arterial whole blood may be used. The system is suitable for use in hospitals, doctor's offices and nursing homes. The HemoCue Glucose 201+ system is for *In Vitro* Diagnostic use only. The results can be printed or downloaded to a PC.

⚠ All system components are designed and manufactured to provide maximum safety. Any other use of the system may impair the safety.

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Components
Bestandteile
Composants
Onderdelen

GB

1. The HemoCue Glucose 201+ Analyzer.*
2. AC adapter.**
3. 4 type AA or R6 batteries, 1.5 V.***
4. A vial of HemoCue Glucose 201 Microcuvettes. Individually packaged HemoCue Glucose 201 Microcuvettes.***
5. HemoCue Glucose 201+ Operating manual and Quick Reference Guide.
6. HemoCue Cleaner.

The HemoCue Glucose 201+ and its AC adapter are delivered in a carton. Open the carton on a stable surface and lift out the analyzer and accessories. Consult local environmental authorities for proper disposal of batteries.

- * Do not open the cover of the analyzer.
- ** Δ Only use adapters as listed under Specifications.
- *** Not included.

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
GB

- 1a. If mains power is available, connect the adapter to the socket on the back.
- 1b. If mains power is not available, insert the 4 type AA or R6 batteries, 1.5 V, into the battery compartment.
- 1c. If a battery symbol appears on the display the batteries are running low. The analyzer will continue to give accurate results but the batteries should be replaced as soon as possible.
2. The analyzer can be connected to a PC (for further information, see **Set-up** Audio signal, time and date or separate manual for PC connection) or directly to a printer (see **Set-up** Printer function).
3. Pull the cuvette holder out to its loading position.
4. Press and hold the left button until the display is activated (all symbols appear on the display).
5. The display shows the version number of the programme, after which it will say "GLU" and "GLU". During this time the analyzer will automatically verify the performance of the optronic unit.
6. After 10 seconds the display will show three flashing dashes and the HemoCue symbol. This indicates that the HemoCue Glucose 201+ is ready for use.

To turn off the analyzer press and hold the left button until the display reads OFF and becomes blank.

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
GB

The time and date function within the set-up menu must be activated before data can be transferred to a PC.

1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
- 3ab. Use the right button to scroll until the display shows an audio symbol in the upper right corner. The signal can be activated or deactivated by pressing the left button.
4. Continue to scroll using the right button until the display shows characters for time, date and year. The hour figure will be flashing.
5. Use the right button to change between hours, minutes, year, day and month. To change a flashing figure hold the left button down for quick advancement or advance slowly by pressing and quickly releasing the button.
6. When all the settings are completed, hold the right button down for approximately 5 seconds. The analyzer will then automatically return to the measuring position.

The analyzer is supplied with default time and date settings. These settings can only be changed by using the PC program (for further information, see separate manual for PC connection).

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
Meting Capillair bloed

GB

Note: In cases of severe hypotension or peripheral circulatory failure, glucose measurement from capillary samples may be misleading. In such circumstances it is recommended that the glucose level should be measured using venous or arterial whole blood¹.

1. After start-up, the cuvette holder should be in its loading position. The display will show three flashing dashes and the HemoCue symbol.
2. Make sure the patient's hand is warm and relaxed. Use only the middle or ring finger for sampling. Avoid fingers with rings on.
3. Clean with disinfectant and allow to dry.
4. Using your thumb, lightly press the finger from the top of the knuckle towards the tip. This stimulates the blood flow towards the sampling point.

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

GB

5. For best blood flow and least pain, sample at the side of the fingertip, not in the centre.
6. Whilst lightly pressing towards the fingertip, prick the finger using a lancet.
7. Wipe away the first 2 or 3 drops of blood.
8. Re-apply light pressure towards the fingertip until another drop of blood appears.

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

GB

9. When the blood drop is large enough, fill the cuvette in one continuous process. Do NOT refill!
10. Wipe off excess blood on the outside of the cuvette tip. Make sure that no blood is drawn out of the cuvette during this procedure.
11. Look for air bubbles in the filled cuvette. If present, take a new sample. Small bubbles around the edge can be ignored.

Note: If a second sample is to be taken from the same fingerstick, it is important that this should be done after the first sample has been read. Wipe away the remains of the first drop of blood and take a second sample from a new drop of blood.

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
GB

12. Place the filled cuvette in the cuvette holder. This should be performed within 40 seconds after filling the cuvette!
13. Push the cuvette holder to its measuring position.
14. During the measurement "⌘" will be shown on the display.
15. After 40–240 seconds the glucose value of the sample is displayed. The result will remain on the display as long as the cuvette holder is in the measuring position. When operating on battery power the analyzer will automatically turn off after approximately 5 minutes if its not used.*
16. Although the reagents are present in the cuvette in extremely low quantities, consult local environmental authorities for proper disposal. Always handle blood specimens with care, as they might be infectious.

* All results are stored in the analyzer and can be viewed on the display using the scroll function (see **Set-up** Memory function – scroll) or by using the PC program (for further information, see separate manual for PC connection).

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Measuring Venous or arterial blood

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


④

GB

1. Mix the blood well before performing the measurement.
- 2ab. Place a drop of blood on a hydrophobic surface, e.g. a plastic film, using a pipette or a DIFF-Safe™.
3. Fill the cuvette in one continuous process. Do NOT refill! Wipe off excess blood on the outside of the cuvette tip. Make sure that no blood is drawn out of the cuvette during this procedure.
4. Look for air bubbles in the filled cuvette. If present, take a new sample. Small bubbles around the edge can be ignored.

Perform the analysis as per steps 12–16 in **Measuring** Capillary blood.

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
Set-up QC-test

GB

The HemoCue® system can be controlled daily with whole blood, haemolysate, serum or waterbased controls. To obtain reliable results only use controls recommended by HemoCue.

1. Press both buttons at the same time.
2. The display shows a flashing QC-symbol.
3. Select QC-test by pressing the left button.
4. The analyzer automatically returns to its measuring position and the QC-symbol appears on the display.
5. Fill the cuvette with control solution recommended by HemoCue and carry out the measurement in accordance with steps 12–16 in **Measuring** Capillary blood.
6. Once the measurement has been made the analyzer automatically returns to its measuring position and the QC-symbol disappears from the display. Patient sampling can begin.

If more QC-tests are to be performed, proceed in the same way following steps 1–6. To deactivate the QC-test, follow steps 1–2 and scroll using the right button until another set-up activity is shown on the display. Hold the right button down for approximately 5 seconds. The analyzer automatically returns to its measuring position and the QC-symbol is no longer visible.

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
Set-up Geheugenfunctie – Scrollen

GB

The analyzer automatically stores up to 600 results. When the memory is full the analyzer will automatically overwrite the oldest result. To view results stored in the memory use the following procedure.

- 1ab. The cuvette holder may be in the loading position or in the measuring position. If in the measuring position the most recent measurement is displayed.
2. Press the left or right button to scroll backwards or forwards between the results. The stored values are now shown on the display. When the button is held down the user can see the order of the results (the most recent measurement is always number 1) and when the button is released the result is shown.
- 3ab. If no buttons are pressed while in the loading position the analyzer automatically resets after 5 seconds and new measurements can be made. If in the measuring position pull out the cuvette holder and wait for the HemoCue symbol after which a new measurement can be made.

Note: The analyzer will not respond to other commands while a measurement is being made. If the error code E35 is shown while scrolling through the results, the measurement has not been stored in the memory, see the **Troubleshooting guide**.

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
Set-up Resultaten wissen

GB

1. Press both buttons at the same time.
2. The display now shows a flashing QC-symbol.
3. Scroll using the right button until the display shows a flashing waste bin in the left corner.
- 4ab. To delete the most recent result press and hold the left button. A warning triangle together with a "1" appears on the display for approximately 3 seconds. The most recent result has now been deleted and the display shows a "0" to confirm that this has been done. The most recent result can only be deleted immediately after the measurement has been made.
- 5ab. To delete all stored results hold both buttons down at the same time for 5 seconds. The display shows the waste bin, the warning triangle and the total number of results stored. Continue to hold both buttons down until the number changes to "000" and the warning triangle disappears. All stored results have now been deleted.

Note: Deleted results cannot be recovered. If results are to be saved they must be printed out or downloaded to a PC before starting the deleting process.

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

Configurer Mode impression

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GB

The analyzer can be on or off when connected to the printer. Only the most recent result can be printed out directly after the measurement.

1. Connect the cable* to the analyzer and printer*.
2. Perform the analysis by following the steps 12–16 in **Measuring** Capillary blood.
3. When the glucose result is shown on the display, the printer will automatically print the result, date and time.

* Not included.

Note: Only use ASCII printers recommended by HemoCue.

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Maintenance Maintenance quotidienne

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GB

The cuvette holder should be cleaned after each day of use.

1. Pull the cuvette holder out to the loading position.
2. While pressing the catch, carefully rotate the cuvette holder sideways in open position as far as possible to the left.
3. Remove the cuvette holder from the analyzer, it will come off the stainless steel pin it rotates on.
4. Clean the cuvette holder with alcohol (20-70%) or mild detergent

Wait 15 minutes before replacing the cuvette holder and using the analyzer. Make sure the cuvette holder is dry before inserting.

The optical parts should be cleaned when directed to do so in the Troubleshooting Guide, see **Maintenance** Optical parts.

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

Maintenance Optical parts

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GB

Dirty optical parts may cause an error code. Follow step 1-3 under Daily maintenance and then clean optical parts as follows. Note! Make sure that the HemoCue® Cleaner reaches both upper and lower cover glasses, see pictures 5a and 5c.

5. With the cuvette holder removed from the analyzer push the Cleaner into the opening of the optic unit, as far in as possible. Move from side to side 5-10 times, see picture 5b and thereafter push in and pull out the Cleaner 5-10 times, cleaning the cover glasses, placed to the left, see picture 5c. If the Cleaner is stained, repeat with a new Cleaner.
6. Wait 15 minutes before replacing the cuvette holder and using the analyzer. Make sure the cuvette holder is dry before inserting.

The cover may be cleaned with alcohol (20-70%) or mild detergent.

Note! As an alternative to the HemoCue® Cleaner, a lint free cotton swab, non-pretreated, moistened with alcohol (20-70% without additive) or water may also be used. If a cotton swab is used make sure it is not too wet and not too dry. Use a dry swab to wipe away excess liquid in the optic house after cleaning with a moistened swab. To avoid scratches on the cover glasses, only the cotton part of the swab should come in contact with the cover glasses.

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GB Troubleshooting Guide

If you are unable to resolve the problem by following this Troubleshooting Guide, please contact your local HemoCue distributor or HemoCue AB. Note! Do not open the cover of the analyzer.

Symptom	Explanation	Action
The analyzer shows an error code.	May be a temporary fault.	Turn off the analyzer and turn it on again after 30 seconds. Take a new microcuvette and repeat the measurement. If the problem continues, see specific error code below.
E00	No stable endpoint is found within the time range. 1. The cuvette is faulty. 2. The circuit board is out of order.	1a. Check the expiry date for the microcuvettes. 1b. Take a new microcuvette and repeat the measurement. 2. The analyzer needs service. Contact your distributor.
E01–E05	1. Fault or dirt in the optronic or electronic units.	1a. Turn off the analyzer and clean the optronic unit using the HemoCue Cleaner. 1b. The analyzer needs service. Contact your distributor.
E06	1. Unstable blank value. The analyzer might be cold.	1. Turn off the analyzer and allow it to reach room temperature. If the problem continues, the analyzer needs service. Contact your distributor.
E07	1. The battery power is too low.	1a. The batteries need to be replaced. Turn off the analyzer and replace the batteries, 4 type R6 or AA, 1.5 V. 1b. Use the mains power supply.
E08	The absorbance is too high. 1. Light blocking item in the cuvette holder.	1a. Check that the analyzer and microcuvettes are being used according to the HemoCue Glucose 201+ operating manual and instructions for use. 1b. The analyzer needs service. Contact your distributor.
E09–E30	1. Fault or dirt in the optronic or electronic units.	1a. Turn off the analyzer and clean the optronic unit using the HemoCue Cleaner. 1b. The analyzer needs service. Contact your distributor.
E35	1. The current measurement has been incorrectly stored and will remain as error code E35 until the entire memory has been deleted.	1a. No action is normally required. 1b. If storage problems continue, the analyzer needs service. Contact your distributor.

Symptom	Explanation	Action
<p>E70–E71</p> <p>This function is only available from serial number 1035 142 xxx</p>	<p>Erroneous kinetics.</p> <ol style="list-style-type: none"> 1. The cuvette is faulty. 2. The sample might be grossly lipemic. 	<ol style="list-style-type: none"> 1-2. Check that the system is used according to the instruction for use. <p>Fill a NEW cuvette and place it immediately in the analyzer.</p> <p>If the error code appears again, use a suitable laboratory reference method to analyze the specimen.</p>
<p>ECU</p> <p>This function is only available from serial number 0835 142 xxx</p>	<ol style="list-style-type: none"> 1. The cuvette is empty or not filled with enough sample. Please note, the ECU function captures primarily empty cuvettes, not under-filled cuvettes. 2. No chemical reaction is identified in the blood filled cuvette. 	<ol style="list-style-type: none"> 1-2. Fill a NEW cuvette and place it immediately in the analyzer. <p>If the ECU function appears again, contact your local distributor or HemoCue AB. Use another HemoCue analyzer or a suitable laboratory reference method to analyze the specimen.</p>
<p>HHH</p>	<ol style="list-style-type: none"> 1. Measured value exceeds 22.2 mmol/L (400 mg/dL). 	<ol style="list-style-type: none"> 1. The measuring range may be extended to 44.4 mmol/L (800 mg/dL) by dilution, see Specifications Results.
<p>No characters on the display.</p>	<ol style="list-style-type: none"> 1. The analyzer is not receiving power. 2. If on battery power, the batteries need to be replaced. 3. The display is out of order. 	<ol style="list-style-type: none"> 1a. Check that the AC adapter is connected to the mains power supply. 1b. Check that the AC adapter is securely connected to the analyzer. 1c. Check that the cable is not damaged. 2. Turn off the analyzer and replace the batteries, 4 type R6 or AA, 1.5 V. 3. The analyzer needs service. Contact your distributor.
<p>The display gives erroneous characters.</p>	<ol style="list-style-type: none"> 1. The display is out of order. 2. The microprocessor is out of order. 	<ol style="list-style-type: none"> 1. The analyzer needs service. Contact your distributor. 2. The analyzer needs service. Contact your distributor.

Symptom	Explanation	Action
The display shows "🔋".	<ol style="list-style-type: none"> 1. The batteries need to be replaced. 2. If on mains power, the AC adapter or the circuit board is out of order. 	<ol style="list-style-type: none"> 1. Turn off the analyzer and replace the batteries, 4 type R6 or AA, 1.5 V. 2a. Check that the AC adapter is properly connected and working. 2b. The analyzer needs service. Contact your distributor.
The display does not switch from "🔋" and "GLU" to three flashing dashes and "👤" (ready for measuring).	<ol style="list-style-type: none"> 1. The magnet in the cuvette holder may be missing. 2. The magnetic sensor is out of order. 	<ol style="list-style-type: none"> 1. The analyzer needs service. Contact your distributor. 2. The analyzer needs service. Contact your distributor.
Measurements on control materials are out of range – either too HIGH or too LOW.	<ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. The control has not been mixed properly and/or is not at room temperature. 4. Air bubbles in the microcuvette. 5. The optronic unit is dirty. 6. The control is not suitable for use with the HemoCue Glucose 201+ system. 7. The calibration of the analyzer has been changed. 8. The controls are beyond their expiry dates or have been improperly stored. 	<ol style="list-style-type: none"> 1. Check the expiry date and the storage conditions of the microcuvettes. 2. Remeasure the sample with a new microcuvette. 3. Make sure that the control is mixed properly and at room temperature. 4. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 5. Clean the optronic unit using the HemoCue Cleaner. 6. Only use controls intended for the HemoCue Glucose 201+ system recommended by HemoCue. 7. The analyzer needs service. Contact your distributor. 8. Check the expiry date and the storage conditions of the control. Take a new microcuvette and repeat the measurement.
Measurements on patient samples are higher or lower than anticipated.	<ol style="list-style-type: none"> 1. The microcuvettes are beyond their expiry date, damaged or have been improperly stored. 2. The optical eye of the microcuvette is contaminated. 3. Air bubbles in the microcuvette. 4. The optronic unit is dirty. 5. The calibration of the analyzer has been changed. 	<ol style="list-style-type: none"> 1. Check the expiry date and the storage conditions of the microcuvettes. 2. Remeasure the sample with a new microcuvette. 3. Check the microcuvette for air bubbles. Remeasure the sample with a new microcuvette. 4. Clean the optronic unit using the HemoCue Cleaner. 5. The analyzer needs service. Contact your distributor.

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Specifications

General

The HemoCue® Glucose 201+ is a system for the determination of the total amount of glucose in whole blood. The system consists of a specially designed analyzer with specially designed cuvettes containing dried reagents. The cuvette serves as a pipette, reaction vessel and as a measuring cuvette. No dilution is required. The glucose measurement takes place in the analyzer, which follows the progress of the reaction and presents the result only when the end point of the reaction has been reached. The system is factory calibrated and needs no further calibration. The HemoCue glucose reference system is traceable to an Isotope Dilution Gas Chromatography - Mass Spectrometry (ID GC-MS) method.

Intended purpose

Quantitative determination of glucose in whole blood using a specially designed analyzer, the HemoCue® Glucose 201+. The quantitative determination of the instant blood glucose concentration in circulation supplements the clinical evidence in the diagnosis and treatment of the diabetic patient as well as in the monitoring of neonatal blood glucose levels. To establish HemoCue Glucose reference values and an intervention level, neonatal blood samples should be evaluated against a suitable laboratory method, taking into consideration the difference between whole blood and plasma reference values. HemoCue Glucose 201 Microcuvettes are for In Vitro Diagnostic use only.

The HemoCue Glucose 201+ Analyzer is only to be used with HemoCue Glucose 201 Microcuvettes. For professional use only.

IVD Medical Device Directive

The HemoCue® Glucose 201+ complies with the IVD Medical Device Directive 98/79/EC and carries the CE mark.

Theory

The chemistry method utilised by the HemoCue® Glucose 201

Microcuvette is a modified glucose dehydrogenase method described by Banauch et al.² A chromogen compound is added to the reagents according to the principle outlined by Bergmeyer³ with saponin used for haemolysing the erythrocytes. The absorbance is measured at two wavelengths (667 and 840 nm) to compensate for turbidity.

Reagents

Saponin, NAD, MTT and NaF.

Enzyme mix: Glucose Dehydrogenase, Diaphorase and Mutarotase.

Specimen collection and preparation

Capillary, venous or arterial blood may be used. The anticoagulants EDTA, sodium heparin and lithium heparin with and without gel and the glycolysis inhibitors sodium fluoride, sodium oxalate and potassium oxalate shall be used. Glycolysis is a major concern in all glucose measurements. To minimise the effect of glycolysis, measure the blood sample as soon as possible after taking the sample. Samples of blood collected in containers with recommended anticoagulants should be analysed within 30 minutes. Mix all samples thoroughly by inverting the sample-container at least 10 times before measurement. Note! Lithium heparin with gel can not be used after centrifugation.

Storage and environmental requirements for HemoCue® Glucose 201 Microcuvettes

Use the HemoCue® Glucose 201 Microcuvettes prior to their expiry date. The expiry date is printed on each package.

Storage for cuvettes kept in a vial

Store unopened HemoCue® Glucose 201 Microcuvettes below 8 °C (46 °F) (incl. storage in a freezer). Note, microcuvettes stored in a freezer must be allowed to reach room temperature (approx. 30 minutes) before analysis. Cuvettes kept in an opened vial are stable for 30 days when stored in a refrigerator at 2–8 °C (35 – 46 °F).

After breaking the seal, the vial can be stored at room temperature for up to 3 days. Close the lid immediately after cuvettes are removed from the vial.

Storage for individually packaged cuvettes

Store unopened HemoCue® Glucose 201 Microcuvettes below 8 °C (46 °F) (incl. storage in a freezer). Note, microcuvettes stored in a freezer must be allowed to reach room temperature (approx. 30 minutes) before analysis. The individually packaged cuvettes can be stored at room temperature for up to 3 days.

HemoCue® Glucose 201+ Analyzer

The analyzer can be stored at 0–50 °C (32–122 °F). Operating temperature is 15–30 °C (59–86 °F). Allow the analyzer to reach ambient temperature before use. The analyzer should not be operated at high (i.e. > 90 % non-condensing) humidity.

Quality Control

The HemoCue® Glucose 201+ Analyzer has an internal quality control, the "selftest". Every time the analyzer is turned on it will automatically verify the performance of the optronic unit of the analyzer. This test is performed every second hour if the analyzer is left turned on. If quality control checks are required for regulatory reasons they should be performed using liquid controls recommended by HemoCue.

Results

The measuring range is 0–22.2 mmol/L (0–400 mg/dL). Results above 22.2 mmol/L (400 mg/dL) will be displayed as HHH. Any results with HemoCue® Glucose Systems suggesting clinical intervention in the hyperglycemic range on pre-term neonates (<37 weeks), should be verified against a suitable laboratory reference method. The measuring range may be extended to 44.4 mmol/L (800 mg/dL) by dilution with saline 1+1. Mix the

sample thoroughly before dilution and measurement. The result on the display should be multiplied by a factor of 2.

Note: Dilution may reduce accuracy. For specific performance characteristics, see the HemoCue Glucose 201 Microcuvettes instructions for use.

Limitations

HemoCue® Glucose 201 Microcuvettes are for In Vitro Diagnostic use only. The HemoCue Glucose 201+ Analyzer is only to be used together with HemoCue Glucose 201 Microcuvettes. For further limitations of the procedure, see the HemoCue Glucose 201 Microcuvettes instructions for use.

Expected values

Fasting glucose values, whole blood, adults 3.5– 5.3 mmol/L (65–95 mg/dL).⁴ For diagnosis of diabetes mellitus, follow local recommendations or use the following value according to WHO: Fasting whole blood glucose, capillary or venous \geq 6.1 mmol/L (\geq 110 mg/dL).⁵

Technical Specifications

Dimensions: 85 × 160 × 43 mm

Weight: 350 g (batteries included)

4 batteries type AA or R6, 1.5 V

AC adapter: CE marked

Only use adapters as listed in the table at the end of

Specifications.

Pollution degree: 2

Overvoltage category: II

The instrument is tested according to IEC 61010-1, First edition, 1990 with Amendment 1, 1992 and Amendment 2, 1995 (EN 61010-1: 1993 with A2: 1995), EN 60601-1-2 and complies with the IVD Medical Device Directive 98/79/EC.

Warranty

The analyzer carries a 24-month warranty from the day of receipt. After the warranty period service/repair is carried out at fixed prices.

Service and Disposal

The analyzer should be cleaned as recommended under Maintenance prior to service or disposal. Consult your local environmental authorities for proper disposal.

Spare parts – Accessories

The following spare parts and accessories are available:

AC adapter

Battery Lid

Cuvette holder

PC program

Cables for PC and printer connection

Liquid controls

HemoCue Cleaners

Hard Carrying Case

HemoCue Lancets

Recommended separation distance between Portable and mobile RF communications equipment and HemoCue® Glucose 201* Analyzer

The HemoCue systems are intended for use in an electromagnetic environment in which radiated RF disturbances are controlled.

The customer or the user of HemoCue systems can help prevent electromagnetic interference by maintaining a minimum distance between portable and RF communications equipment (transmitters) and HemoCue systems as recommended below, according to the maximum output power of the communications equipment.

Rated maximum output power of transmitter (W)	Separation distance according to frequency of transmitter (m)		
	150 kHz to 80 MHz $d=1.2\sqrt{P}$	80 MHz to 800 MHz $d=1.2\sqrt{P}$	800 MHz to 2.5 GHz $d=2.3\sqrt{P}$
0.01	0.12	0.12	0.23
0.1	0.38	0.38	0.73
1	1.2	1.2	2.3
10	3.8	3.8	7.3
100	12	12	23

For transmitters rated at maximum output power not listed above, the recommended separation distances (d) in meters (m) can be estimated using the equation applicable to the frequency of the transmitter, where (P) is the maximum output power rating of the transmitter in watts (W) according to transmitter manufacturer.

NOTE 1 At 80 MHz and 800 MHz, the separation distance for the higher frequency range applies.

NOTE 2 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

Guidance and manufacturer's declaration – Electromagnetic immunity

The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or user of the HemoCue systems should assure that it is used in such an environment.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
Electrostatic discharge (ESD) IEC 61000-4-2	±6 kV contact ±8 kV air	±6 kV contact ±8 kV air	Floors should be wood, concrete or ceramic tile. If floors are covered with synthetic material, the relative humidity should be at least 30 %.
Electrical fast transient/ burst IEC 61000-4-4	±2 kV for power supply lines ±1 kV for input/output lines	2 kV for power supply lines ±1 kV for input/output lines	Mains power quality should be that of a typical commercial or hospital environment.
Surge IEC 61000-4-5	±1 kV differential mode	±1 kV differential mode	Mains power quality should be that of a typical commercial or hospital environment.
Voltage dips, short interruptions and voltage variations on power supply lines IEC 61000-4-11	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds For explanation of U see NOTE 1	<5 %U (>95 % dip in U) for a 0.5 cycle 40 %U (60 % dip in U) for 5 cycles 70 %U (30 % dip in U) for 25 cycles <5 %U (>95 % dip in U) for 5 seconds	Mains power quality should be that of a typical commercial or hospital environment. If the user of the HemoCue systems requires continued operation during power mains interruptions, it is recommended that the HemoCue systems be powered from an uninterruptible power supply or a battery.

Immunity test	IEC 60601 test level	Compliance level	Electromagnetic environment – guidance
<p data-bbox="302 256 431 303">Conducted RF IEC 61000-4-6</p> <p data-bbox="302 338 431 385">Radiated RF IEC 61000-4-3</p>	<p data-bbox="540 256 708 303">3 Vrms 150 kHz to 80 MHz</p> <p data-bbox="540 338 708 385">3V/m 80 MHz to 2.5 GHz</p> <p data-bbox="540 421 751 438">See NOTE 2 and NOTE 3</p>	<p data-bbox="777 256 846 274">3 Vrms</p> <p data-bbox="777 338 837 356">3 V/m</p>	<p data-bbox="950 256 1338 415">Portable and mobile RF communications equipment should be used no closer to any part of the HemoCue systems, including cables, than the recommended separation distance calculated from the equation applicable to the frequency of the transmitter</p> <p data-bbox="950 421 1227 438">Recommended separation distance</p> <p data-bbox="950 473 1024 491">$d=1.2\sqrt{P}$</p> <p data-bbox="950 509 1190 526">$d=1.2\sqrt{P}$ 80 MHz to 800 MHz</p> <p data-bbox="950 544 1190 562">$d=2.3\sqrt{P}$ 800 MHz to 2.5 GHz</p> <p data-bbox="950 597 1325 697">Where (P) is the maximum output power rating of the transmitter in watts (W) according to the transmitter manufacturer and (d) is the recommended separation distance in meters (m).</p> <p data-bbox="950 732 1338 832">Field strengths from fixed RF transmitters, as determined by an electromagnetic site survey (a), should be less than the compliance level in each frequency range (b).</p> <p data-bbox="950 867 1325 909">Interference may occur in the vicinity of equipment marked with the following symbol.</p> <div data-bbox="1109 938 1175 1003" style="text-align: center;">  </div>

Technical specifications (EMC-RF)

Use only cables with the following specification:

USB shielded maximum 2 m

Serial shielded maximum 1.5 m

Guidance and manufacturer's declaration – electromagnetic emissions		
The HemoCue systems are intended for use in the electromagnetic environment specified below. The customer or the user of the HemoCue systems should assure that it is used in such an environment.		
Emission test	Compliance	Electromagnetic environment – guidance
RF emissions	Group 1	The HemoCue systems uses RF energy only for its internal function. Therefore, its RF emissions are very low and not likely to cause any interference in nearby electronic equipment.
RF emissions	Class B	The HemoCue systems are suitable for use in all establishments, including domestic establishments and those directly connected to the public low-voltage power supply network that supplies buildings used for domestic purposes.
Harmonic emissions IEC 61000-3-2	Class A	
Voltage fluctuations/flicker emissions IEC 61000-3-3	Complies	

NOTE 1 It is the manufacturer's responsibility to provide equipment electromagnetic compatibility information to the customer or user.

NOTE 2 It is the user's responsibility to ensure that a compatible electromagnetic environment for the equipment can be maintained in order that the device will perform as intended.

NOTE 1 U is the a.c. mains voltage prior to application of the test level.

NOTE 2 At 80MHz and 800MHz, the separation distance for the higher frequency range applies.

NOTE 3 These guidelines may not apply in all situations. Electromagnetic propagation is affected by absorption and reflection from structures, objects and people.

- a) Field strengths from fixed transmitters, such as base stations for radio (cellular/cordless) telephones and land mobile radios, amateur radio, AM and FM radio broadcast and TV broadcast cannot be predicted theoretically with accuracy. To assess the electromagnetic environment due to fixed RF transmitters, an electromagnetic site survey should be considered. If the measured field strength in the location in which the HemoCue systems are used exceeds the applicable RF compliance level above, the HemoCue systems should be observed to verify normal operation. If abnormal performance is observed, additional measures may be necessary, such as re-orienting or relocating the systems.
- b) Over the frequency range 150 KHz to 80 Mhz, field strength should be less than 3 V/m.

Symbols used

Caution

Consult instructions for use

CE mark

Class II equipment

Only valid within the European Community.
Indicates separate collection for waste of electrical and electronic equipment.

Temperature limitation

Efficiency Level

To maintain safety use only adapter marked HCA01

DC Inlet

Relative humidity, non-condensing

References

1. Atkin et al, Annals of Internal Medicine, 1991, 114;12, 1020–1024.
2. Banauch et al, Z. Klin. Chem. u. Klin. Biochem, 1975;13:101–107.
3. Bergmeyer, Methods of Enzymatic Analysis 1, 1974; Chemie Publishers, Weinheim.
4. Tietz Fundamentals of Clinical Chemistry: Carbohydrates, Fifth Edition, Editor, W.B. Saunders Company, 2001, 427–461.
5. Definition and Diagnosis of Diabetes Mellitus and Intermediate Hyperglycaemia. Report of a WHO/IDF Consultation. Geneva: WHO 2006.

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HemoCue AB, Box 1204, SE-262 23 Ängelholm, Sweden. Phone +46 77 570 02 10. Fax +46 77 570 02 12.